

Gingival Health

in vivo study

Gingival health of a periodontal maintenance gingival health assessment after use of Sonicare Elite

Donly K. University of Texas. Data on file

Objective	To demonstrate the ability of Sonicare Elite to reduce gingivitis.
Methodology	One hundred subjects participated in this parallel, single-blind study. Subjects were randomly assigned to treatment groups within strata defined by high or low plaque scores at screening. Subjects were evaluated for gingivitis at baseline and after two, four and twelve weeks of product use. Gingivitis was assessed with the Löe and Silness Gingival Index.
Results	Subjects using Sonicare Elite showed a significant reduction in gingival index scores in diseased sites after 12 weeks compared to baseline ($p < 0.0001$).
Conclusion	Sonicare Elite reduces gingivitis.

